

Rapid Differentiation of Isobaric Drugs Using a Novel Direct Sample Analysis Source, CID, and Accurate Mass TOF Mass Spectrometry

Robert J. Seward, Bonnie Marmor, and Andrew Tyler

1 Introduction

Identification of isobaric drugs is an increasing challenge for analytical research laboratories. These isomers cannot be distinguished by molecular weight alone. Liquid or gas chromatographic methods can often separate structural isomers, but sample throughput may be reduced. Collisionally induced dissociation (CID) in the mass spectrometer has the potential to distinguish structural isomers without the need for chromatography. Here, we examine the potential to differentiate isobaric drugs using direct sample ionization combined with accurate mass TOF MS analysis.



Figure 1. AxION 2 TOF MS fitted with a Direct Sample Analysis (DSA) Source.

- Analysis was performed on a PerkinElmer AxION 2 TOF MS system fitted with a Direct Sample Analysis (DSA) source operated in positive mode.
- Samples were acquired using AxION DSA Controller and data was processed with AxION Solo software using a strong signal setting of 5000 counts and a cutoff of 1% strong signal except where noted. Isotope search window ± 0.003 Da. Monoisotopic weight 7.
- Analytes were obtained from Cerilliant and Sigma-Aldrich.
- Analytes were prepared in methanol to 1 $\mu\text{g/mL}$ and 5 μL was spotted on the stainless steel sample mesh.
- 5 replicates were acquired per sample, 1 second acquisition per sample, 10 spectra per second acquisition rate.

3 Results




Figure 2. Example analysis of synthetic cannabinoid mixture containing an isobaric pair: JWH-019 and JWH-122 which were distinguished by unique fragment ions. Identified ions are highlighted in color.




Figure 3. Analysis of the isobaric pair cocaine and scopolamine. (a) Cocaine only. (b) Scopolamine only. (c) 1:1 cocaine and scopolamine mix. Identified ions are highlighted in color.

