The Uncertainty of Reference Standards

A Guide to Understanding Factors Impacting Uncertainty, Uncertainty Calculations and Vendor Certifications

Authors: Ning Chang PhD, Isil Dilek PhD, Kevin Gates, Huahua Jian PhD, Sherri Pogue, Uma Sreenivasan PhD

© 2009 Cerilliant Corporation



Introduction

Analytical results are only as good as the calibrators used in the analysis.

Certified ampouled solution reference standards are widely used in the forensic/toxicology, clinical/diagnostic, and pharmaceutical industries to support forensic investigations, therapeutic monitoring, and clinical decisions.

Accuracy is dependent on: robustness of the analytical method, preparation of samples and standards, & the quality, purity and accuracy of the reference.



Uncertainty

Uncertainty is "a parameter associated with the result of a measurement that characterizes the dispersion of the values that could reasonably be attributed to the measurand." *

*International Vocabulary of Basic and General Terms in Metrology

For solution reference standards, whether purchased or prepared at the bench, <u>the measurand in question is the</u> <u>concentration</u> of the solution (specified in units of mass per volume).

Understanding uncertainty variables is important for laboratories seeking to comply with ISO/IEC 17025 requirements and for those preparing reference solutions from neat materials at the bench.



Uncertainty is an estimation

Factors Impacting Concentration

- Neat material purity & homogeneity
- Solution preparation techniques/methods
- Weighing equipment & procedures
- Pipettes and flasks
- Solution storage, handling, & homogeneity



Cerilliant's Approach to Reference Solution Preparation

- Complete characterization of neat materials
- Qualified and controlled weighing processes
- Gravimetric addition of solvent
- Sealed into ampoules under argon
- Analytical verification of concentration, purity, & homogeneity
- Shelf life determined through real-time stability studies



Methodology for Uncertainty

Kragten Spreadsheets and differential perturbation

W. Guthrie, T. Vetter. "Hands-on Workshop on Evaluating Uncertainties for Chemical Analysis" Gaithersburg, MD: National Institute of Science and Technology, PITTCON 2007.

Cerilliant's Measurement Equation

$$C = \frac{(m_{v+a} - m_v)dp}{(m_{f+s} - m_f)} \pm U$$

Where: C = Concentration of solution (mass/volume) $m_{v+a} = mass of analyte + vial$

 $m_v = mass of empty vial$

 $m_{f+s} = mass of flask + solvent$

- $m_f = mass of empty flask$
- d = density of solution
- U = the assigned expanded measurement uncertainty



Factors Impacting Uncertainty of Solution Standard Preparations



Uncertainty Associated with Purity of Neat Material

- Is the neat material certification adequate for reference standard use? – Includes Purity and Residuals?
- The Cerilliant Purity Factor mass balance equation often referred to as "potency" on vendor COA's – is used to calculate the amount of material needed to achieve accurate concentration of the solution standard.

 $PurityFactor = \left[[100 - (wt\%OVI) - (wt\%H_2O) - (wt\%ROI)] * \frac{ChromPurity}{100} \right] \pm U$

wt%OVI: the weight percentage of residual solvents present in the neat material.

wt%H2O: the weight percentage of water present in the neat material.

wt%ROI: the weight percentage of inorganic content in the neat material.

ChromPurity: based on the chromatographic purity of the specified primary purity method, either GC or HPLC.

U = the assigned expanded measurement uncertainty





Key Elements of Uncertainty Neat Material Purity Factor



Purity Factor Considerations

- Appropriate methods for purity determination
- Quantitation of residual impurities
 - Water
 - Solvent
 - Inorganics
- Use of a mass balance equation
- Stability purity <u>and</u> residuals



Factors Impacting Uncertainty of Solution Standard Preparations Every step of the process has Neat Material uncertainty and must be evaluated **Purity Factor** $u_{pf} = 0.292\%$ Chromatographic Purity Solvent Addition Solution Density **Residual Water Analysis Residual Solvent Analysis** Mass measurement Inorganic Content Analysis Temperature Uncertainty Instrument Tolerances Of Solution Weighing Techniques Concentration Preparation Steps **Balance Sensitivity & Linearity** Balance Selection, Qualification – Minimum Weights Mass Measurement 11

Uncertainty of the Mass Measurement *u*_m

Factors impacting mass measurement uncertainty

- Mass measurement uncertainty included the following components:
 - *u*_{sens} Uncertainty due to the balance's sensitivity tolerance. The
 sensitivity tolerance includes the uncertainty of the balance's built-in
 reference weight used for the internal calibrations.
 - u_{lin} Uncertainty due to non-linearity of the characteristic curve.
 - s_p Repeatability includes effects from readability, drift, static, ambient drafts, thermal drafts, vibration, gross/net weight, eccentric loading, temperature stability, EMI/RFI, weighing procedure, installation, tare container geometry, adsorption/absorption, and balance settings.



Repeatability is specific to YOUR environment, YOUR balances, and YOUR weighing technique

Mass Measurement Considerations

Appropriate balance selection and qualification are critical to ensuring accuracy of the solution standard and can have a significant impact on the overall uncertainty

 Improper balance selection can lead to high levels of uncertainty

Importance of Balance Selection and Mass Uncertainty				
Mass Uncertainty				
	5-place Balance	4-place Balance		
1 mg	8.0%	45.0%		
10 mg	0.80%	4.5%		
100 mg	0.080%	0.45%		
1000 mg	0.0080%	0.045%		

 Minimum weighings established to achieve USP specified minimum relative error of NMT 0.1%.

Cerilliant Minimum Weighing Requirements

Balance	7-place	6-place	5-place	4-place
Balance Resolution	0.0001 mg	0.001 mg	0.01 mg	0.1 mg
Minimum Weighing	1 mg	3 mg	20 mg	125 mg

Mass Measurement Considerations

Weighing Technique can significantly influence uncertainty

- Appropriate assignment of uncertainty of solution standard preparation must consider weighing technique in addition to balance selection and qualification.
- Accuracy of weighing can be influenced by:
 - tongs vs. gloved hands
 - balance equilibration time
 - sample and solvent temperature
 - ambient temperature
 - vibrations
 - movement of air

For Example:

Cerilliant studies indicate that when gloved hands are used as opposed to tongs for handling sample vials, uncertainty of mass measurement increased approximately 10 fold.

• Air currents, drafts around the balance, and additional vibrational forces on the pan can significantly affect balance repeatability.



Uncertainty of Mass Measurement u_m

<u>Repeatability</u> experiments are used to determine preparer specific uncertainty for weighing operations

Repeatability Experiments for Cerilliant Preparations					
		Using Differe	nt Balances		
Balance	XP6400	XP123OS	XP205	XP56	UMX2
Duluille	1 Place	3 Place	5 Place	6 Place	7 Place
Process Scale	1-10 Liters	100-250 mL	25 mg-10 g	5 mg	l mg
Approx. Gross Mass	1 kg	200 g	2.1 g	2.005 g	41 mg
Tare Container	none	none	2 mL glass vial	2 mL glass vial	aluminum micro weigh pan
Ref./Net Mass (g)	1000	200	0.1	0.005	0.001
Mean	1000.000	200.0005	0.099995	0.0050023	0.00100193
Standard Deviation	0.032	0.0012	0.000024	0.0000018	0.00000035
%RSD (s _p)	0.00324%	0.00060%	0.02350%	0.03500%	0.0352248%



Uncertainty of the Mass Measurement v_m

• The largest standard deviation observed in repeatability tests per process was used to calculate the uncertainty value for weighing operations

 $(u_m = 0.035\%)$

$$u_m = \sqrt{s_p^2 + u_{lin}^2 + u_{sens}^2}$$

• The term for uncertainty from the weighing operation is **applied to each weighing operation** in the measurement equation for solution standard preparation, mass measurement of the analyte and mass measurement of the solvent.

$$C = \frac{(m_{v+a} - m_v)dp}{(m_{f+s} - m_f)} \pm U$$



Solvent Addition Temperature vs. Density

Change in density with temperature can affect volumetric preparation of a solution but can be controlled by gravimetric addition of solvent



Bench preparation of sample and reference on different days may create variability due to density change 0.57% difference in concentration when prepared volumetrically at 20° vs. 25°C





Source: Handbook of Thermophysical and Thermochemical Data, CRC Press

Addition of Solvent / Density Uncertainty U_d

- Gravimetric addition controls for variation of solvent density (and therefore volume) with temperature
 - traceability to SI units of mass.
 - Solvent mass required is calculated from the density
 - Solvent is added by weight
 - Solution mass is converted back to volume by dividing the mass by the density.

(For low concentration solutions (< 2 mg/mL), the solution density is approximated using the density of the pure solvent)

- Density is measured using a Mettler Toledo Densito 30PX density meter which has a resolution of 0.0001 g/mL.
 - The uncertainty component for density was estimated based on instrument tolerances

$$u_d = \frac{0.001}{\sqrt{3}} = 0.000577 \text{ g/mL}$$



Calculation of the Combined Standard Uncertainty (*u_c*) and Expanded Uncertainty (*U*) of the Solution Standard Concentration

Range of Variables Modeled				
Variable/Result Name	Lowest Modeled Value	Highest Modeled Value		
Analyte mass	l mg	5 g		
Solvent mass	60 g	1.2 kg		
Adjustment factor	1	1.01		
Solution density	0.6 g/mL	1.2 g/mL		
Batch volume	100mL	1000mL		
Concentration	1 μg/mL	5 mg/mL		

Inputs to the measurement equation for concentration could take on a wide range of values depending on batch volume, target concentration, solution density, and purity of analyte. Input values were, therefore, varied to provide models which yielded uncertainty values for at least 99% of gravimetrically prepared solution.

Over the range of models tested the relative expanded uncertainties varied little, reinforcing the value and importance of process controls employed.



Cerilliant Model

Uncertainty of Solution Standard Preparations



Considerations

4 Points of Consideration

- 1. Thorough characterization of the neat material is essential to determine an accurate mass balance purity factor.
 - Characterization of the neat material should be appropriate for use
 - Should include both purity and residuals
 - Uncertainty contributions from neat material certification are significant
- 2. Actual practice in the laboratory can vastly influence uncertainty related to weighing operations.
 - Balance manufacturer specifications are insufficient
 - Repeatability assessment must be included
- 3. Gravimetric preparations provide greater control and traceability for standard preparation.
 - Volumetric dilutions should account for errors arising temperature/density effects and user error associated with visual read lines.
- 4. Variations or changes to any component of the solution standard preparation process can impact uncertainty and requires reassessment of uncertainty values.



Considerations

Evaluating Certificates of Analysis

- How was the uncertainty determined? What quality systems were used (such as ISO/IEC 17025 or ISO guide 34)?
- What does the uncertainty value cover?
- Is the uncertainty reported as an expanded uncertainty with a coverage factor? Are confidence intervals provided?
- Is the neat material traceability and test data provided?
- Is purity of the neat material considered in the uncertainty of the standard preparation? Was the purity method appropriate for the compound, sufficiently robust and repeatable?
- What components are included in the Purity Factor assessment? Were residuals considered? What methods were used to determine these values?
- Are environmental conditions such as temperature or density considered in the uncertainty statement?
- Are balance and volumetric tolerances included those of the manufacturer alone or experimentally verified for the manufacturing process?



If the standard was diluted from a stock solution, does the uncertainty include uncertainty of the secondary dilution?



A-064 FN022009-02 Revision 0 Page 1 of 5

Certificate of Analysis

ISO/IEC 17025

ISO GUIDE 34

ISO 9001:2000

Acetaminophen 4-Acetamidophenol

Catalog Number:	A-064	H
Solution Lot:	FN022009-02	
Expiration Date:	February 2014	
Solvent:	Methanol	
Volume per Ampule:	Not less than 1 mL	HU I
Storage:	Protect from air and light, refrigerate or freeze.	
Intended Use:	For laboratory use only. Not suitable for human or appinal consumption.	

· Expiration Date has been established through real time stability studies.

 Ampules are overfilled to ensure a minimum 1 mL volume fill. We advise laboratories to use measured volumes of this standard solution before diluting to the desired concentration.

Comp	onent	Chromatographic Parity	Concentration		
Acetar	ninophen	99.9%	1.000 ± 0.006 mg/mL€		
Uncertainty of the concentration is expressed as an expanded uncertainty in accordance with ISO 17025 and Guide 34 at the approximate					
95% cor	95% confidence interval using a coverage factor of k = 2 and has been calculated by statistical analysis of our production system and				
incorpor	incorporates uncertainty of the purity factor, material density, and mass.				
C	Consideration is a smooth different consideration of deal and so with a first structure of an ideal is smoothing				

Solution Standard Verification and Homogeneity

Standard		Verified Concentration (mg/mL)		%RSD - Homogeneity	
Solution	Lot Number	Actual Results	Acceptance Crtieria	Actual Results	Acceptance Crtieria
New Lot	FN022009-02	0.998	± 3%	0.8	≤ 3%
Previous Lot	FN092707-01	0.996	± 3%	0.4	≤ 3%

Concentration is verified through multiple analyses and is calculated as the average of multiple analyses compared to an independently prepared calibration curve

 Homogeneity of the New Lot is ensured through rigorous production process controls statistically analyzed to evaluate risk and verified by analysis. The % RSD of samples pulled from across the lot demonstrate homogeneity of the New Lot.

The % RSD of the Previous Lot represents variability of the analysis

Traceability

+ Gravimetrically prepared using qualified balances calibrated semi-annually by Mettler Toledo using NIST traceable weights. Calibration verification performed weekly and prior to each use utilizing NIST traceable weights. Each balance has been assigned a minimum weighing by Mettler Toledo taking into consideration the balance and installed environmental conditions to ensure weighing complies with USP tolerances of no more than 0.1% relative error.

Concentration is verified against an independently prepared 4-point calibration curve gravimetrically prepared using balances calibrated to NIST

Cerilliant certifies that this standard meets the specifications stated in this certificate and warrants this product to meet the stated acceptance criteria through the expiration/retest date.

ACLASS

Cerilliant Corporation

Lara Sparks, Quality Assurance Director

March 18, 2009 Date

Solution Purity is verified chromatographically post ampouling to ensure no degradation or contamination

Concentration & Uncertainty of the gravimetric preparation expressed as:

1.000 + 0.0006 mg/mL

Description of Cerilliant's Uncertainty value & confidence interval:

"Uncertainty of the concentration is expressed as an expanded uncertainty in accordance with ISO 17025 and Guide 34 at the approximate 95% confidence interval using a coverage factor of k = 2 and has been calculated by statistical analysis of our production system and incorporates uncertainty of the purity factor, material density, and mass".

Analytical Verification of Concentration & Homogeneity

Gravimetrically prepared concentration is verified analytically. Acceptance criteria incorporates variability of the analysis. Homogeneity is verified analytically by analyzing ampoules pulled from across the lot.

Traceability Statement describing traceability to SI units

"Gravimetrically prepared using qualified balances calibrated semiannually by Mettler Toledo using NIST traceable weights. Calibration verification performed weekly and prior to each use utilizing NIST traceable weights. Each balance has been assigned a minimum weighing by Mettler Toledo taking into consideration the balance and installed environmental conditions to ensure weighing complies with USP tolerances of no more than 0.1% relative error".



A-064 FN022009-02 Revision 0 Page 2 of 5

Calibration Curve

Linear Regression

4

0.999

Calibration Curve:

Number of Points:

Linearity (r):

Corilliant
Analytical Reference Standards
,

Neat Material Characterization Summary and Purity Factor Assignment

Shows detail of all analyses performed and results - allows confirmation of Purity Factor variables considered

811 Paloma Drive, Suite A, Round Rock, TX 78665 Cerilliant Corporation

Standard Solution Assay Parameters

HPLC/UV Analysis Method: Column: Luna 5µ C18, 4.6 x 250 mm Mobile Phase: Methanol::Water (50::50) Flow Rate: 1.0 mL/min Wavelength: 246 nm

Neat Material Data

Compound Name: Compound Lot:	Acetaminophen PN061507-01	Chemical Formula: CAS Number: Molecular Weight:	C ₈ H ₉ NO ₂ 103-90-2 151.16
	Neat Material Cha.	racterization Summary	
Analytical Test		Method	Results
Primary Chromatographic Purity b	y HPLC/PDA Analysis	SP10-0102	99.9%
Secondary Chromatographic Purity	y by GC/FID Analysis	SP10-0101	99.9%
Identity by GC/MS Analysis		SP10-0105	Consistent with Structure
Identity by ¹ H-NMR Analysis		SP10-0116	Consistent with Structure
Residual Solvent Analysis by GC/	FID Headspace	AM1087	0.00%
Residual Water Analysis by Karl F	ischer Coulometry	SP10-0103	< 0.2%
Inorganic Content by Microash Ar	alysis	Outsourced	< 0.1%
Purity Factor			99.90%

· Primary purity is calculated as the average of two independently performed analyses utilizing two different methods. Acceptance criteria requires the purity values to be within 0.5% of each other.

· The primary chromatographic purity value is used to calculate the Purity Factor.

· A secondary chromatographic purity method is utilized as a control.

Purity Factor = [(100 - wt% residual solvent - wt% residual water - wt% residual inorganics) x Chromatographic Purity/100].

Purity factor does not include adjustment for chiral and/or isotopic purity.

Acknowledgements and References

Acknowledgements:

• Ning Chang PhD, Isil Dilek PhD, Kevin Gates, Huahua Jian PhD, Sherri Pogue, Kristine Waddell

<u>References:</u>

- ISO/IEC10725:2005, 2nd ed., "General Requirements for the competence of testing and calibration laboratories"
- ISO GUIDE 34:2000, 2nd ed., "General Requirements for the competence of reference material producers" and ISO Guide 35:2006, 3rd ed., "Reference materials-General and statistical principles for certification"
- International Organization for Standardization, "International Vocabulary of Basic and General Terms in Metrology", International Organization for Standardization, 1993
- EURACHEM/CITAC Guide 2nd ed., "Quantifying Uncertainty in Analytical Measurement", EURACHEM/CITAC, 2000
- a) W. Guthrie, T. Vetter. "Hands-on Workshop on Evaluating Uncertainties for Chemical Analysis" Gaithersburg, MD: National Institute of Science and Technology, PITTCON 2007. b)T. Vetter. Quantifying Measurement Uncertainty in Analytical Chemistry – A Simplified Practical Approach, Measurement Science Conference 2001, Session V-B, Anaheim, CA, January 19, 2001 web link: <u>http://www.cstl.nist.gov/acd/839.03/Uncertainty.pdf</u>
- J. Kragten. Calculating Standard Deviations and Confidence Intervals with a Universally Applicable Spreadsheet Technique. *The Analyst* 119: 2161-2165 (1994)
- Arthur Reichmuth, "The Uncertainty of Weighing Data Obtained with Electronic Analytical Balances," Mettler Toledo 2004







science, smarter.

Cerilliant Quality

ISO GUIDE 34

ISO/IEC 17025

ISO 9001:2000

GMP/GLP