

# Evaluation of LC-MS/MS Scrambling Ratios for Deuterium-Labeled Vitamin D Metabolites, Steroids and Other Compounds of Clinical Significance

Authors: [Joshua Cooper](#), Huahua Jian, Derrell Johnson, Isil Dilek, and Uma Sreenivasan  
Cerilliant Corporation, 811 Paloma Drive, Suite A, Round Rock, TX 78665

## Abstract

**Introduction and Objective:** A significant clinical challenge with LC-MS/MS is the potential for matrix effects that cause interferences or impact ionization efficiency. Stable isotope-labeled internal standards are frequently used to compensate for matrix effects and to increase the accuracy of quantitation. The use of a labeled internal standard that co-elutes with the drug being monitored can potentially offset patient specific matrix effects (co-eluting concomitant medication, etc.) that may occur at the retention time of the analyte of interest. Complications in the use of deuterium-labeled internal standards can arise from hydrogen-deuterium scrambling in the collision cell at the selected transitions or in the ion source. In this study, we examined deuterium labeled 25-Hydroxyvitamin D, testosterone, and other compounds of clinical significance by LC-MS/MS at multiple transitions. We investigated reproducibility of the scrambling ratio and influences on scrambling of different LC-MS systems (tandem quadrupole vs. quadrupole time-of-flight), matrix selection, concentration, and deuterium placement in the internal standard.

## Methods and Procedures

LCMS System 1:

Instrument: Waters Alliance UPLC-Xevo G2 Q-ToF  
Column: Waters Acquity UPLC, BEH C18, 1.7 $\mu$ m, 2.1 x 50mm

25-Hydroxyvitamin D Analysis Conditions:

UPLC Conditions: 0.4mL/min, gradient, 0.1:99.9 to 99.9:0.1 (0.1% formic acid in acetonitrile:0.1% formic acid in water)  
MS Conditions: ESI+, Cone 25V, Capillary 2.5kV, CE 20

Testosterone Analysis Conditions:

UPLC Conditions: 0.4mL/min, isocratic, 30:70 (0.1% formic acid in acetonitrile:0.1% formic acid in water)  
MS Conditions: ESI+, Cone 30V, Capillary 3.0kV, CE 18

LCMS System 2:

Instrument: Agilent 1100 HPLC-6410 triple quad  
Column: Phenomenex Kinetex, C18, 3 $\mu$ m, 2.1 x 50mm

25-Hydroxyvitamin D Analysis Conditions:

HPLC Conditions: 0.4mL/min, isocratic, 80:20 (0.1% formic acid in methanol:0.1% formic acid in water)  
MS Conditions: ESI+, Fragmentor 110V, Capillary 4.0kV, CE 5

Testosterone Analysis Conditions:

UPLC Conditions: 0.4mL/min, isocratic, 30:70 (0.1% formic acid in acetonitrile:0.1% formic acid in water)  
MS Conditions: ESI+, Fragmentor 50V, Capillary 4.0kV, CE 10


Solution Standards Used:

25-Hydroxyvitamin D3, Cat# H-083  
25-Hydroxyvitamin D3-d<sub>6</sub>, Cat# H-074  
25-Hydroxyvitamin D2, Cat# H-073  
Testosterone, Cat# T-037  
Testosterone-d<sub>3</sub>, Cat# T-046  
Testosterone-<sup>13</sup>C<sub>3</sub>, Cat# T-037  
Progesterone-d<sub>6</sub>, Cat# P-070  
Pregabalin-d<sub>6</sub>, Cat# P-072


Serum Extraction:  
200 $\mu$ L of sample in serum + 200 $\mu$ L of methanol, vortexed to mix.  
Added 1mL of heptane, vortexed for 30sec,  
Centrifuged for 4min at 3000rpm  
900 $\mu$ L of top layer dried under nitrogen  
Reconstituted in 100 $\mu$ L of ethanol

## Comparisons of 25-Hydroxyvitamin D2 and D3 Deuterium Scrambling


25-Hydroxyvitamin D2




25-Hydroxyvitamin D2-d<sub>3</sub>




25-Hydroxyvitamin D2-d<sub>6</sub>



25-Hydroxyvitamin D3




25-Hydroxyvitamin D3-d<sub>6</sub>




Labeled 25-Hydroxyvitamin D2 and D3 Scrambling in Serum

| Compound             | Label | System  | Concentration $\mu$ g/mL | Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ |
|----------------------|-------|---------|--------------------------|-----------------------|-----------------------|------------------------------|
| 25-Hydroxyvitamin D2 | $d_3$ | Xevo G2 | 2                        | 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 28.6                         |
|                      |       |         | 0.2                      | 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 35.4                         |
|                      |       |         | 5                        | 416 $\rightarrow$ 397 | 416 $\rightarrow$ 398 | 2.8                          |
|                      |       |         |                          | 416 $\rightarrow$ 379 | 416 $\rightarrow$ 380 | 19.7                         |
|                      |       |         |                          | 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 30.4                         |
|                      |       |         |                          | 416 $\rightarrow$ 397 | 416 $\rightarrow$ 398 | 2.8                          |
|                      | $d_6$ | 6410    | 50                       | 416 $\rightarrow$ 379 | 416 $\rightarrow$ 380 | 20                           |
|                      |       |         |                          | 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 30.5                         |
|                      |       |         |                          | 419 $\rightarrow$ 400 | 419 $\rightarrow$ 401 | 2                            |
|                      |       |         |                          | 419 $\rightarrow$ 382 | 419 $\rightarrow$ 383 | 8.8                          |
|                      |       |         | 5                        | 401 $\rightarrow$ 382 | 401 $\rightarrow$ 383 | 5.9                          |
|                      |       |         |                          | 419 $\rightarrow$ 400 | 419 $\rightarrow$ 401 | 2                            |
| 25-Hydroxyvitamin D3 | $d_6$ | 6410    | 50                       | 419 $\rightarrow$ 382 | 419 $\rightarrow$ 383 | 9                            |
|                      |       |         |                          | 401 $\rightarrow$ 382 | 401 $\rightarrow$ 383 | 5.4                          |
|                      |       |         |                          | 407 $\rightarrow$ 388 | 407 $\rightarrow$ 389 | 4                            |
|                      |       |         |                          | 407 $\rightarrow$ 370 | 407 $\rightarrow$ 371 | 18.8                         |
|                      |       |         | 2.5                      | 389 $\rightarrow$ 370 | 389 $\rightarrow$ 371 | 9.2                          |
|                      |       |         |                          |                       |                       |                              |

Vitamin D in Serum on 6410



Vitamin D in EtOH Scrambling on Xevo G2



Transitions Comparisons for Native and Labeled 25-Hydroxyvitamin D2 and D3 in EtOH on 6410

| Parent $\rightarrow$ Water loss |        |                          |                       |                       |                              |
|---------------------------------|--------|--------------------------|-----------------------|-----------------------|------------------------------|
| Compound                        | Label  | Concentration $\mu$ g/mL | Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ |
| 25-Hydroxyvitamin D2            | $d_3$  | 100                      | 416 $\rightarrow$ 397 | 416 $\rightarrow$ 398 | 2.9                          |
|                                 | $d_6$  | 100                      | 419 $\rightarrow$ 400 | 419 $\rightarrow$ 401 | 2                            |
|                                 | native | 50                       | 413 $\rightarrow$ 394 | 413 $\rightarrow$ 395 | 0.5                          |
| 25-Hydroxyvitamin D3            | $d_6$  | 50                       | 407 $\rightarrow$ 388 | 407 $\rightarrow$ 389 | 4                            |
|                                 | native | 100                      | 401 $\rightarrow$ 382 | 401 $\rightarrow$ 383 | 0.5                          |

Parent  $\rightarrow$  2 Water losses

| Compound             | Label  | Concentration $\mu$ g/mL | Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ |
|----------------------|--------|--------------------------|-----------------------|-----------------------|------------------------------|
| 25-Hydroxyvitamin D2 | $d_3$  | 100                      | 416 $\rightarrow$ 379 | 416 $\rightarrow$ 380 | 19.5                         |
|                      | $d_6$  | 100                      | 419 $\rightarrow$ 382 | 419 $\rightarrow$ 383 | 8.9                          |
|                      | native | 50                       | 413 $\rightarrow$ 376 | 413 $\rightarrow$ 377 | 0.5                          |
| 25-Hydroxyvitamin D3 | $d_6$  | 50                       | 407 $\rightarrow$ 370 | 407 $\rightarrow$ 371 | 18.9                         |
|                      | native | 100                      | 401 $\rightarrow$ 364 | 401 $\rightarrow$ 365 | 0.3                          |

Water Loss  $\rightarrow$  2 Water losses

| Compound             | Label  | Concentration $\mu$ g/mL | Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ |
|----------------------|--------|--------------------------|-----------------------|-----------------------|------------------------------|
| 25-Hydroxyvitamin D2 | $d_3$  | 100                      | 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 30.4                         |
|                      | $d_6$  | 100                      | 401 $\rightarrow$ 382 | 401 $\rightarrow$ 383 | 5.4                          |
|                      | native | 50                       | 398 $\rightarrow$ 376 | 398 $\rightarrow$ 377 | 0.4                          |
| 25-Hydroxyvitamin D3 | $d_6$  | 50                       | 389 $\rightarrow$ 370 | 389 $\rightarrow$ 371 | 11.2                         |
|                      | native | 100                      | 383 $\rightarrow$ 364 | 383 $\rightarrow$ 365 | 0.3                          |

Notes: 25-Hydroxy D2-D6 water loss $\rightarrow$ 2 water loss has same transition as 25-Hydroxyvitamin D3 parent $\rightarrow$ water loss. Can be problem if compounds are not well resolved chromatographically

Selection of Transitions Greatly Impacts Observed Scrambling


5 $\mu$ g/mL Infusion at 20 $\mu$ L/min of  $d_3$  labeled 25-Hydroxyvitamin D2 on Xevo G2

| Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ |
|-----------------------|-----------------------|------------------------------|
| 416 $\rightarrow$ 397 | 416 $\rightarrow$ 398 | 2.2                          |
| 416 $\rightarrow$ 379 | 416 $\rightarrow$ 380 | 16.9                         |
| 398 $\rightarrow$ 379 | 398 $\rightarrow$ 380 | 30.9                         |


Note: Under optimized UPLC-Q-ToF conditions only water loss MS ions were detected. MS ion ratios changed for 25-Hydroxyvitamin D when combined with mobile phase. Could detect ions without water loss when infusing.

## Investigation of Testosterone Scrambling


Testosterone




Testosterone-<sup>13</sup>C<sub>3</sub>



Testosterone-d<sub>3</sub>



Testosterone Chromatograms on 6410




Testosterone Scrambling Comparison


| Label  | Method   | Instrument | Concentration $\mu$ g/mL | Transitions $D_{n-1}$ or $^{13}C_{n-1}$ | Transitions $D_n$ or $^{13}C_n$ | *Scrambling % $D_{n-1} / D_n$ |
|--------|----------|------------|--------------------------|-----------------------------------------|---------------------------------|-------------------------------|
| $d_3$  | Infusion | Q-ToF      | 10                       | 292 $\rightarrow$ 255                   | 292 $\rightarrow$ 256           | 31.9                          |
|        |          |            | 100                      |                                         |                                 | 36.5                          |
|        |          |            | 10                       |                                         |                                 | 35.7                          |
|        | LC       | 6410       | 100                      |                                         |                                 | 37.7                          |
|        |          |            | 10                       |                                         |                                 | 36.3                          |
|        |          |            | 100                      |                                         |                                 | 0.1                           |
| native |          |            | 100                      | 289 $\rightarrow$ 252                   | 289 $\rightarrow$ 253           | 0.0                           |

Major transitions are:  
Native: 289 $\rightarrow$ 97 & 289 $\rightarrow$ 109  
Testosterone-d<sub>3</sub>: 292 $\rightarrow$ 97 & 292 $\rightarrow$ 109  
Testosterone-<sup>13</sup>C<sub>3</sub>: 292 $\rightarrow$ 100 & 292 $\rightarrow$ 112  
No scrambling at major transitions

Testosterone Scrambling at m/z 253



Testosterone Scrambling at 97 and 109




Testosterone  $d_{n-2} / d_n$  Scrambling

| Label | Method   | Instrument | Concentration $\mu$ g/mL | Transition $d_{n-2}$  | Transition $d_n$      | Scrambling % $d_{n-2} / d_n$ |
|-------|----------|------------|--------------------------|-----------------------|-----------------------|------------------------------|
| $d_3$ | Infusion | Q-ToF      | 10                       | 292 $\rightarrow$ 254 | 292 $\rightarrow$ 256 | 2.6                          |
| $d_3$ | LC       | Q-ToF      | 100                      | 292 $\rightarrow$ 254 | 292 $\rightarrow$ 256 | 3.6                          |
| $d_3$ | LC       | Q-ToF      | 10                       | 292 $\rightarrow$ 254 | 292 $\rightarrow$ 256 | <LOD                         |


## Scrambling for other clinical compounds

Xevo G2 Scrambling Infusion Experiments


| Compound     | Label | Transition $d_{n-1}$  | Transition $d_n$      | Scrambling % $d_{n-1} / d_n$ | Transition $d_{n-1}$ |
|--------------|-------|-----------------------|-----------------------|------------------------------|----------------------|
| Progesterone | $d_6$ | 324 $\rightarrow$ 305 | 324 $\rightarrow$ 306 | 20                           | 19                   |
|              |       | 324 $\rightarrow$ 287 | 324 $\rightarrow$ 288 | 77                           | 19                   |
|              |       | 324 $\rightarrow$ 112 | 324 $\rightarrow$ 113 | 0                            | 19                   |
|              |       | 324 $\rightarrow$ 99  | 324 $\rightarrow$ 100 | 0                            | 19                   |
| Pregabalin   | $d_6$ | 166 $\rightarrow$ 147 | 166 $\rightarrow$ 148 | 0                            | 25                   |
|              |       | 166 $\rightarrow$ 129 | 166 $\rightarrow$ 130 | 0                            | 25                   |
|              |       | 166 $\rightarrow$ 102 | 166 $\rightarrow$ 103 | 12                           | 25                   |
|              |       | 166 $\rightarrow$ 88  | 166 $\rightarrow$ 89  | 40                           | 25                   |



Progesterone-d<sub>6</sub>



Pregabalin-d<sub>6</sub>



## CONCLUSIONS

- Scrambling was observed on both the Agilent 6410 triple quadrupole and the Waters Xevo G2 Q-ToF, and in some cases was very pronounced.
- For a specific transition, scrambling ratios were consistent between solvent and serum. No matrix effects on scrambling.
- Direct infusion can provide rapid and accurate determination of scrambling ratios. Infusion and chromatographic injection results were consistent.
- It may be advisable to investigate at higher concentrations than normally analyzed to ensure that instrument sensitivity does not impact accuracy of scrambling determination.
- Awareness of potential scrambling is important for proper internal standard selection. Scrambling may be mitigated or eliminated by altering instrument conditions and transition selection.
- Deuterium-labeled internal standards are a viable option for LC-MS/MS analysis with selection of the appropriate transition. Deuterated standards can be more cost effective than <sup>13</sup>C labeled internal standards, more widely available and with lower cost per test. <sup>13</sup>C labeled internal standards are most effective when deuterium scrambling issues can not be resolved.